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Setting the Matrix Printer Epson LX-300+

The matrix printer Epson 300+ provides a serial interface RS 232 for the connection to BINDER units with RD2 controller. The printer will print out the temperature data in intervals that can be set at the temperature chamber controller.

In order to adapt the printer to the data transfer protocol of the RD2 controller, a few parameter changes have to be carried out at the printer as described in this manual.

The parameter changes follow a menu which will be printed out automatically during the setting procedure. Each menu point is encoded by a combination of LED lights,

e.g.,

The icons symbolize the LED status:

- LED off
- LED on
- LED flashes

The first two symbols refer to the font LEDs, the third one is the PAUSE LED. The LED status therefore indicates, which is the actually entered menu point or parameter.


Figure 1: Printer display

In a first step, the language is chosen. Following this step, a menu overview with information about the individual setting options will be printed in the selected language.

The necessary parameters of the printer setting are completely described in this manual, so it is not necessary to rely on the printout.

1. Preparation

Charge the printer with paper according to the printer's operating manual instructions. Switch on the printer at the main switch while pressing down the TEAR OFF key.

The first printout will be printed, containing multi-language instructions how to select the language of further printouts, as well as an explanation of the LED symbols used. Right after this printout, you have to carry out the language settings.

2. Language setting

You have the choice of the following languages:

- ■ □ English
- □ ■ Français
- □ □ Deutsch
- ■ ■ Italiano
- ■ □ Español
- □ ■ Português

Select the desired language by generating the following combination of LED signals via the LF/FF key:

- for language setting German: **off – on – on**
- for language setting English **off – off – on**
- for language setting French **off – on – off**

(proceed accordingly to select other languages)

Confirm via the TEAR OFF key.

Following your choice of language (if needed, after releasing the PAUSE button), an overview of the printer setting options will be printed in the selected language.

** Menu **

- ■ □ Page length for tractor
- □ ■ Skip over perforation
- □ □ Auto tear off
- ■ ■ Auto line feed
- ■ □ Print direction
- □ ■ Software
- ■ ☼ 0 slash
- ☼ ■ High speed draft
- ☼ ☼ I/F mode
- ☼ ■ ■ Auto I/F wait time
- ☼ ■ ☼ Baud rate
- ☼ ☼ ■ Parity
- ■ ☼ Data length
- ☼ ■ Parallel I/F bidirectional mode
- ☼ ☼ Packet mode
- □ ☼ Character table
- ☼ □ International character set for Italic table
- ☼ ■ □ Manual feed wait time
- ☼ □ ■ Buzzer
- ☼ □ □ Auto CR (IBM 2380 Plus)
- ☼ ☼ □ IBM character table

Pressing upon the TEAR OFF key will confirm your intention to go through the setting procedure.

A few pages menu overview with information about the individual setting options will now be printed.

3. Interface settings (I/F mode)

To enter the sub-menu „I/F mode“, generate the following combination of LED signals via the TEAR OFF key:

off – flash – flash

and confirm via the LF/FF key.

The menu offers the following selection:

```
** I/F mode **  
■ ■ □   Auto  
■ □ ■   Parallel  
■ □ □   Serial
```

Select the Setting „Serial“ by generating the following combination of LED signals via the LF/FF key:

off – on – on

and confirm via the TEAR OFF key.

4. Setting the Baud rate

To enter the sub-menu „Baud rate“, generate the following combination of LED signals via the TEAR OFF key:

flash – off – flash

and confirm via the LF/FF key.

The menu offers the following selection:

```
** Baud rate **  
■ ■ □   19200 BPS  
■ □ ■   9600 BPS  
■ □ □   4800 BPS  
□ ■ ■   2400 BPS  
□ ■ □   1200 BPS  
□ □ ■   600 BPS  
■ ■ ☼   300 BPS
```

Select the Setting „9600 BPS“ by generating the following combination of LED signals via the LF/FF key:

off – on – off

and confirm via the TEAR OFF key.

In order to save the changes, switch off the printer at its main switch.

5. Settings at the RD2 controller

To specify the print interval, proceed as follows:

Switch off the temperature controller by pressing key .

Press Mode-key for 10 about seconds, until the gradient of the set-point or “r OFF“ appears on the display. Press the Mode-key another two times.

The setting of the print interval P appears on the display. Default setting is P 60, which means one printout every 60 minutes. Set the printing interval to the desired value (e.g. P1, that is one printout every minute) by keys . After about 10s the controller returns to normal display.

Start up the controller by pressing key .


Figure 2: Controller display RD2

6. Putting the printer into operation

Switch on the printer at the main switch, insert paper and press the form feed key LF/FF, until the paper covers the roller.

The font can be selected according to the Epson 300+ operating manual (press down the PAUSE key for 3 sec, then choose the font via the TEAR OFF key). The choice of a large font, e.g., EPSON Sans Serif or EPSON Roman, is recommended.

The actual temperature shown on the display is printed according to the selected print interval, e.g. “Temp.: 100“ on the left paper edge. A few cm beside it is printed the sign * to graphically show the continuation of the temperature. Every 5th printout, the printing interval („Ptime: 001“) appears instead of the numerical output. But the * sign is printed as before in value-depending distance.

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- LED **arrêt**
- LED **marche**
- ☼ LED **clignotant**

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