

Steps to select Sterilizer;

Intended use

- 1. What type of load needs to be processed?**
Load / types / programs / resist hot/ Humidity/ Chemicals
- 2. What type of sterilizer is required?**
Dry hot / EO / Steam / Formaldehyde
- 3. What models are available?**
Manual / Automatic / doors
- 4. Where will the sterilizer be sited?**
Upstairs/ which floor / location in the room Middle / Corner
- 5. What services are requirements?**
External Steam / water / drainage / compressed air / electricity / ventilation
- 6. What size and how many?**
one bigger / two smaller / weekly calculation
- 7. What other equipment will be needed?**
Steam generator / air compressor / De-mineral water
- 8. What are the specifications?**
Which standard / any additional details
- 9. What sort of contract?**
Delivery / installation / inspection / Tester
- 10. Which manufacturer?**
What do you think?
- 11. What happens after delivery?**
Warranty / validation / Service